

# Effectiveness Other Swabs Can't Touch.



## FOCUS ON SANITATION AND DISINFECTION

Now, more than ever, it is crucial to verify the effectiveness of your cleaning and sanitation processes. That means that detection is an even more critical component to healthcare and food safety operations.

### A3 TECHNOLOGY

Officials around the world are emphasizing the need for enhanced disinfection to control the presence of the SARS-CoV-2 virus on surfaces. It is important to know that your heightened processes are being effective and, while Kikkoman A3 Technology (ATP+ADP+AMP) does not directly detect the SARS-CoV-2 virus – *very few technologies can do that* – independent testing has shown that A3 Technology represents a novel ATP test that can detect residues left behind by ineffective cleaning at lower levels than other tests on the market. This makes it an effective means to demonstrate and improve your sanitation and disinfection processes.

**kikkoman®**



## Detection Perfection in 10 Seconds.

### STEP 1 Swab


Swab the sample and insert the swab into main body

### STEP 2 Measure

Shake well and insert into the Lumitester

### STEP 3 Analyze

Measurements (RLU) displayed within 10 seconds



## A New Standard in Clean for Many Sites.

### FOOD SAFETY



### HEALTH CARE



### JANITORIAL SERVICES



### INDUSTRIAL MANUFACTURING



### ■ Rapid Hygiene Monitoring System

Lumitester™ Smart Product code: 61234

Measurement time: 10 seconds

Data output: RLU (Relative Light Unit)

Power: Two AA alkaline or nickel hydride rechargeable batteries

Accessories: Two AA alkaline batteries, cleaning brush, USB cable, strap, Quick Manual

※ Do not use this product for purposes other than hygiene monitoring.  
 ※ It is not to be used for counting general living bacteria or detecting specific pathogens.

### ■ ATP+ADP+AMP Test

LuciPac™ A3 Surface Product code: 60361 100 swabs/kit

LuciPac™ A3 Water Product code: 60365 100 swabs/kit

Storage Condition: 2-8°C (Do not freeze)

25°C: 14 days (Unopened)

30°C: 5 days (Unopened)

Expiry: 15 months after manufacturing date

※ Use LuciPac A3 for Lumitester Smart, PD-30 or PD-20. Do not use it for other models.

FMCG Industry Solutions Pty Ltd  
 Unit 11A, 1 - 3 Endeavour Road,  
 Caringbah NSW 2229  
 Australia

Tel: +61 02 9540 2288  
 Fax: +61 02 9012 0320  
[www.fmcgis.com.au](http://www.fmcgis.com.au)  
[info@fmcgis.com.au](mailto:info@fmcgis.com.au)

