



## Q-Swab™ | Environmental Sample Collection

Hygiena™ Q-Swab is a ready-to-use sample collection and delivery device for environmental surface sampling. Money and time are saved by eliminating pipettes and broth preparation required with other sample collection methods. Snap, squeeze, and swab and your sample is ready.

### Benefits:

- Ready-to-use and convenient
- All-in-one sample collection system
- Self-contained, gamma-irradiated broth in each Q-Swab
- Patented Snap-Valve™ allows quick release of broth
- Rehydrates 1 dried media film perfectly
- Can be used on wet or dry surfaces
- 100% recyclable plastic
- Available broth: Letheen
- 18-month shelf life: Letheen (2 - 25°C)



Catalog No.

Description


Quantity

QS12000

Q-Swab Sample Collection Device, 1.0mL, Letheen

250

## Q-Swab Procedure



### FMCG Industry Solutions Pty Ltd

ABN 22 135 446 007

Factory 11A, 1 - 3 Endeavour Rd, Caringbah NSW 2229

[www.fmcgis.com.au](http://www.fmcgis.com.au)

[sales@fmcgis.com.au](mailto:sales@fmcgis.com.au)

1300 628 104 or 02 9540 2288